



Introducing a  
different type of insulin pump

that's the

twiist™



## The twiist™ Automated Insulin Delivery (AID) System

An **effective** and **more personalized insulin delivery system** for people with type 1 diabetes.<sup>1-6</sup>

# twiist is designed with **confidence in mind**

---

Whether you're currently taking insulin with multiple daily injections or you've been using a pump for a while, it's **time to get to know twiist**. twiist looks and works differently than any insulin pump around.

It offers **unmatched personalization** and breakthrough technology to ensure you, or your loved one, **get the insulin you expect, 24 hours a day**.<sup>1-6</sup>



## Adaptability

### Customize it to your day

With unique features like pre-meal preset, the widest glucose target range, and the ability to go back and change a carb entry or the time you ate, twiist gives you more ways to adapt your pump to your day.<sup>1-6</sup>



## Certainty

### Precision with each dose

twiist is the only pump that measures the actual volume of insulin at every dose to make sure the expected amount of insulin was delivered.



## Effectiveness


### Helps you reach your goals

With twiist, you have the tools you need to help you stay in range and manage your insulin delivery with confidence and control.<sup>7</sup>



# No other device lets you do this much

twiist **adapts to your day** with innovative technology.<sup>1</sup> Take a look at what twiist can do.



## Auto Adjust powered by twiist loop

Every 5 minutes, twiist calculates your 6-hour glucose forecast based on your continuous glucose monitor (CGM) readings, insulin delivery history, carb entries, and personalized settings, and adjusts basal insulin to help keep you in range.




## Pre-meal & activity presets

You can temporarily lower your glucose target before meals and raise it before activity to help prevent you from going too high or low.




## After-meal adjustments

Lets you change the number of carbs or the time you ate, even after a meal. So when plans change, twiist has you covered.



## Absorption rate selector

Lets you specify whether you ate fast-, medium-, or slow-absorbing carbs, and twiist will fine-tune the insulin response to the absorption rate.



## Customizable glucose target(s)

twiist lets you set your target anywhere you need to, from 87-180 mg/dL, the widest glucose target range available.<sup>1-6</sup>



## Wear it how you want

You can easily clip twiist to clothing, carry it in your pocket, or stick it to your skin.\* Move it around depending on your outfit or your plans for the day.



## The only pump with 3 ways to bolus

twiist lets you bolus from the pump or a compatible phone. And it's the only device that also lets you discreetly bolus from a compatible smartwatch.



## Big on convenience, small on size

twiist is the smallest pump to hold 300 units of insulin. And it weighs less than 2 oz.<sup>1-6</sup>


\*Sequel is working on adhesive options for on-body wear. Alternative adhesive options are available through third-party providers.

# Type 1 diabetes management is **always within reach**


## Get to know the **twiist App**

No more juggling multiple devices or apps—the twiist App gives **you a clear view of all your essential data in one place**, including your insulin pump status, CGM readings, and 6-hour glucose forecast.

With the app, you can also customize therapy settings, deliver boluses, log your meals, and activate your pre-meal & activity presets so you can truly personalize how you manage your diabetes with twiist.



You can get both apps from the Apple App Store.




## The **twiist insight™ App** for caregivers

It provides peace of mind for caregivers. The insight App syncs up with the twiist AID System data so you can remotely and securely monitor detailed information, helping you to stay in the know as you work together to manage type 1 diabetes.

Note: The insight App is not intended to be a replacement for the information displayed by the twiist System App and should not inform dosing decisions. It is not intended to replace self-monitoring practices as advised by a physician. Data updates require WiFi or cellular connectivity and their frequency may vary depending on connection strength.

# Two unique technologies. One device.

- **twiist Loop**—predictive and adaptive; adjusts to your day
- **iiSure™**—precise insulin measurement for accurate insulin delivery



## iiSure Technology



## twiist Loop

### iiSure sound wave technology

- Using iiSure technology, twiist is the only pump that measures the **actual volume of insulin down to the microdose** to make sure the expected amount of insulin was delivered<sup>1-6</sup>
- **Detects blockages up to 9x faster**, which can help reduce the risk of unexplained highs<sup>8,9†</sup>

### twiist Loop

- **Based on Tidepool Loop, an intuitive, smart technology** that was originally designed by people with diabetes for people with diabetes
- **Auto Adjust** updates your 6-hour glucose prediction every 5 minutes and adjusts basal insulin all day and night to help prevent highs and lows

<sup>†</sup>When comparing user guides of all commercially available AID systems in the United States as of January 2026 (at 1 U/hr or 2 U/hr).

## Indications for Use

### twiist AID System Indications for Use

The twiist Automated Insulin Delivery (AID) System is for people with type 1 diabetes, 6 years of age and above, available with a prescription, only in the U.S. Visit [twiist.com/safety](https://twiist.com/safety) for important safety information.

### Contraindications

The twiist automated insulin delivery system is not intended for anyone unable or unwilling to:

- Test blood glucose (BG) levels as recommended by their healthcare provider.
- Maintain sufficient diabetes self-care skills.
- See their healthcare provider regularly.
- Demonstrate adequate carbohydrate-counting skills (preferred, not required).

The user must not:

- Have difficulty seeing.
- Have difficulty hearing.
- Have deficient cognitive capabilities.
- Have physical impairments which would make operating the pump or iPhone difficult.


Please visit [twiist.com/safety](https://twiist.com/safety) for important safety information.

Product only available in the U.S.

## References

1. Sequel Med Tech - twiist Automated Insulin Delivery (AID) System Powered by Tidepool User Guide.
2. Insulet - Omnipod 5 Automated Insulin Delivery System User Guide.
3. Medtronic - MiniMed 780G system User Guide.
4. Tandem Diabetes Care - t:slim X2 Insulin Pump User Guide.
5. Tandem Diabetes Care - Tandem Mobi System User Guide.
6. Beta Bionics - iLet Bionic Pancreas System User Guide.
7. Data on file. Sequel Med Tech LLC.
8. Calculation based on a comparison of the following IFUs. The following reported 1 U/hr basal rate: Sequel Med Tech - twiist Automated Insulin Delivery (AID) System Powered by Tidepool User Guide; Insulet - Omnipod 5 Automated Insulin Delivery System User Guide; Medtronic - MiniMed 780G system User Guide; Beta Bionics - iLet Bionic Pancreas System User Guide. The following reported 2 U/hr basal rate: Tandem Diabetes Care - t:slim X2 Insulin Pump User Guide; Tandem Diabetes Care - Tandem Mobi System User Guide.
9. Klonoff DC, et al. Time to Moderate and Severe Hyperglycemia and Ketonemia Following an Insulin Pump Occlusion. *J Diabetes Sci Technol.* 2024;18(6):1472-1479.

<sup>†</sup>The twiist™ Savings Program offers U.S. residents with commercial insurance a \$0 Starter Kit and Refill Kits for \$50 or less (may be higher with high-deductible plans). The maximum benefit is \$200 per Refill Kit. Call 1-877-489-4478 or your insurer to verify coverage.



Sequel is working on adhesive options for on-body wear. Alternative adhesive options are available through third-party providers.

Affordability issues shouldn't get in the way of diabetes management.

that's the **twiist**<sup>™</sup>



Most people start twiist for **\$0 the first month** and no more than **\$50 each month thereafter**.<sup>†</sup> You can even get twiist if you're still under contract with another device.

### Ready to try twiist? Here's what's next:

1. Let your doctor know you're interested in twiist.
2. Give us a call at 1-877-4twiist (1-877-489-4478) Mon-Fri, 9 AM to 9 PM EST for guidance regarding your coverage.
3. Get your prescription twiist Starter Kit from the pharmacy.
4. Arrange your live twiist training session.



twiist works with multiple CGMs, with more compatibility options on the way!

### Find out more

Head to [twiist.com/discover-twiist](https://twiist.com/discover-twiist) for more information



@twiistdiabetes



@twiistdiabetes



@twiistdiabetes



@twiist\_diabetes

"twiist," "iiSure," and the twiist logo are licensed trademarks of Sequel Med Tech, LLC. "Tidepool" and the Tidepool logo are registered trademarks of Tidepool Project, a fully-qualified 501(c)(3) nonprofit organization. All other registered trademarks are the property of their respective owners.

© 2026 Sequel Med Tech, LLC. All rights reserved. MAT-0788 V1.0 02/26